


Labco Screw Cap Vials and Bottles

Dimensions are in mm and include Cap.
Capacities given are brim full in ml.
Diagrams may not be to scale



Vial 4
46 x 15.5mm
4.5ml
Borosilicate Glass
Round Bottom




Vial 1
55 x 15.5mm
5.9ml
Soda Glass
Flat Bottom




Vial 3 Flat
101 x 15.5mm
12.0ml
Soda Glass
Flat Bottom



Vial 3 Borosilicate
101 x 15.5mm
12.0ml
Borosilicate Glass
Round Bottom



Vial RK Borosilicate
104 x 15.5mm
12.0ml
Borosilicate Glass
Round Bottom



Vial CT (Crimp Top)
100 x 15.5mm
11.0ml
Borosilicate Glass
Round Bottom



Vial 2
45 x 12.5mm
3.0ml
Soda Glass
Flat Bottom




Vial 7
54 x 12.5mm
3.7ml
Soda Glass
Flat Bottom



Vial 7R
54 x 12.5mm
3.7ml
Soda Glass
Round Bottom



Vial 6
105 x 12.5mm
7.5ml
Soda Glass
Flat Bottom



Vial 9
46 x 20mm
9.0ml
Soda Glass
Flat Bottom



Vial 8
85 x 28mm
38.0ml
Soda Glass
Flat Bottom